



### HANDLABELER PRINTEX 3 LINES



The handlabeler Printex combined with our documentation labels rolls allows you to print in a fast way all details required from CSSD.

You can print the **LOT number, Sterilizer code, Batch number, Operator, Sterilization Date and expiry Date**. In this way you can have a full track of the sterilization cycle.

**PRODUCT CODE: PC2928EIEEN**

**Labels size (b x h) :** 29 mm x 28 mm

**Label type:** rectangular with single or double adhesive support

**Printing features :** 3 printing lines

**Numbers of digits::** 11 digits per each line.

**Numbers of positions per each digit:** 13

**Digit type:** numbers, letters, currency symbols ,various symbols.

**Digit size (b x h):** 1.56 mm x 3 mm

**Printing size max (b x h):** 23 mm x 3 mm

## Adhesive labels for sterilization department

**Labels roll:** 500 labels /roll

**Material:** ABS, POM, RUBBER, COTTON, INOX STEEL, BRASS, INK

**Dimensions (b x p x h):** 30.00 cm x 16.5 cm x 8.0 cm

**Weight:** 537.00 g

**Dimensions of box (b x p x h):** 27.5 cm x 23.5 cm x 6.5 cm

**Content in 1 box:** Handlabeller with 1 ink roll

### PRODUCT CONFORMITY

Manufactured in Factory with environmental management system certified EN ISO 14001:2004 and Quality Certification EN ISO 9001:2008

### MANUFACTURER

PRINTEX SRL- Colleferro-Roma Italy

