

STERILIZATION REEL



- Medical Paper provides maximum safety.
- All imprint are located outside the packaging area to prevent ink pigment migration to the product.
- Transparent multilayer co-polymer film allows easy identification.
- Reinforced film to avoid tear during opening.
- Superior barrier and quality with medical kraft paper.
- Indicators which are water based non toxic and give accurate result for steam and EO Gas and formaldehyde Sterilization.
- Strong seal-strength, visible sealing lines.
- Clean fibre-free opening.
- All materials comply with international standards.
- Conformity with EN 868-5.